

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N,N'-Hexamethylenebis(methacrylamide) - contains <300 ppm MEHQ as inhibitor

Product Number: **729825**
CAS Number: 16069-15-1



Formula: C₁₄H₂₄N₂O₂
Formula Weight: 252.35 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Light Yellow and Faint Brown to Light Brown
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	62.6 - 70.6 %
Nitrogen	10.4 - 11.8 %
Water (by Karl Fischer)	≤ 6.0 %
Purity (GC)	≥ 94.0 %
MEHQ as Inhibitor	≤ 300 ppm

Specification: PRD.0.ZQ5.10000055837

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.