


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Luperox® LP, Lauroyl peroxide - ≥98%

Product Number: 290785
 CAS Number: 105-74-8
 MDL: MFCD00008964
 Formula: C₂₄H₄₆O₄
 Formula Weight: 398.62 g/mol
 Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Granules
Infrared spectrum	Conforms to Structure
Assay	≥ 98 %
Miscellaneous Assay	≥ 3.93 %
Active Oxygen	
Vendor Information	Confirmed
Product of Arkema Inc.	
Recommended Retest Period	-----
9 Months	

Specification: PRD.6.ZQ5.10000019960

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

