

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

ESCIN

Product Number:

E1378

CAS Number:

6805-41-0

TEST**Specification**

Appearance (Colour)

White to Off White

Appearance (Form)

Powder

Titration with NaOH

95.0 - 105.0 %

(anhydrous)

Purity (TLC)

≥ 95 %

Solubility (Colour)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear to Slightly Hazy

50 mg/ml in Methanol

Water

≤ 5.0 %

Residual Solvents (GLC-HS)

Corresponds to Requirements

< = 2.0 % Ethanol, < = 0.3 % Methanol

Specification: PRE.0.ZQ5.10000021117

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.