

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Heparinase I from Flavobacterium heparinum - Lyophilized powder stabilized with approx. 25% bovine serum albumin,  $\geq 200$  units/mg protein (enzyme + BSA)

Product Number:

**H2519**

CAS Number:

9025-39-2

MDL:

MFCD00131317

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

units/mg protein

 $\geq 200$ 

(Enzyme + BSA)

One unit will form 0.1 micromole of unsaturated uronic acid per hour at pH 7.5 at 25 deg C using Heparin, Sodium as substrate.

Specification: PRD.1.ZQ5.10000022357

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.