

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
FK 269 Co(II) PF₆ salt

Product Number: **805548**

Formula: C₂₂H₁₈CoF₁₂N₁₀P₂
Formula Weight: 771.31 g/mol



TEST

Specification

Appearance (Color)	Faint Orange to Very Dark Orange
Appearance (Form)	Powder or Crystals
Infrared Spectrum	Conforms to Structure
Carbon Content	33.4 - 35.2 %
Nitrogen Content	17.7 - 18.7 %
Product of Supplier	Confirmed
Greatcell Solar®	
Registered Trademark	Note
Greatcell Solar is a registered trademark of Greatcell Solar	
Recommended Retest Period	-----
2 Years	

Specification: PRD.0.ZQ5.10000097671

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.