


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris(triphenylphosphinegold)oxonium tetrafluoroborate

Product Number: 665142
CAS Number: 53317-87-6
MDL: MFCD06658141
Formula: C₅₄H₄₅P₃OAu₃BF₄
Formula Weight: 1,480.56 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Crystal to Powder
Infrared spectrum	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Gravimetric Analysis	38.9 - 40.9 %
%Gold (After Reduction With Sulfite)	
Carbon	42.8 - 44.9 %

Specification: PRD.1.ZQ5.10000045485

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.