


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Itraconazole - ≥98% (TLC)

Product Number: I6657
CAS Number: 84625-61-6
MDL: MFCD00870168
Formula: C₃₅H₃₈Cl₂N₈O₄
Formula Weight: 705.63 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Carbon	58.5 - 60.5 %
Nitrogen	15.5 - 16.2 %
Purity (TLC)	≥ 98 %
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/ml in CHCl ₃	

Specification: PRD.1.ZQ5.10000071108

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.