

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Isopropanol - 70% in H2 O

Product Number: 563935 CAS Number: 67-63-0

MDL: MFCD06412526

Formula Weight: 60.1 g/mol

ОН Н₃С СН₃

TEST

Purity

Specification

Appearance (Color)
Appearance (Form)

%Wt 2-Propanol determined via

difference in water content by KFB.

Expiration Date Period

3 Years

Specification: PRD.2.ZQ5.10000004272

Colorless Liquid 69 - 71 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.