

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-(Trimethylsilyloxy)ethyl methacrylate - contains ≤ 100 ppm BHT as inhibitor, 96%

Product Number:

347485

CAS Number:

17407-09-9

MDL:


MFCD00053869

Formula:

C₉H₁₈O₃Si

Formula Weight:

202.32 g/mol

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

 ≥ 95.5 %

BHT as Inhibitor

 ≤ 100 ppm

Specification: PRD.0.ZQ5.10000020953

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.