

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

**Product Number:****750514**

CAS Number:

308068-56-6



TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Solid
May appear as a thick film powder or sheets	
SWNT	Confirmed
Metal Catalyst	Confirmed
< 1%	
Semi-Conducting Purity	
Specification: Report Result	
Metallic Purity	
Specification: Report Result	

Specification: PRD.0.ZQ5.10000100951

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.