

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Titanium - crystalline, 5-10 mm, ≥99.99% trace metals basis (purity exclusive of Na and K content)

Product Number: 305812
CAS Number: 7440-32-6
MDL: MFCD00011264
Formula: Ti
Formula Weight: 47.87 g/mol

Ti

TEST	Specification
Appearance (Color)	Grey
Appearance (Form)	Crystals
Size (pieces)	≤ 10 mm
Typically 5-10 mm	
ICP Major Analysis	Confirmed
Confirms Titanium Component	
Trace Metal Analysis	≤ 100.0 ppm
Purity	Meets Requirements
≥99.99% Based on Trace Metals Analysis (Exclusive of Na and K Content)	

Specification: PRD.0.ZQ5.10000017253

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.