

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(acrylic acid) – viscosity 350–700 cP , 4.0% mucilage(25 °C)


Product Number:

900338

CAS Number:

9003-01-4

Formula:

(C₃H₄O₂)_n

TEST

Specification

Appearance (Color)

White to Off White

Appearance (Form)

Powder

1H NMR Spectrum

Conforms to Structure

Brookfield Viscosity

350 - 2000 cP

4.0% mucilage, 25°C

Recommended Retest Period

2 Years

Specification: PRD.0.ZQ5.10000093546

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.