


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Benzbromarone – analytical standard

Product Number: B5774
CAS Number: 3562-84-3
MDL: MFCD00078962
Formula: C₁₇H₁₂Br₂O₃
Formula Weight: 424,08 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity)	Clear
100 mg/mL, EtOH	
Carbon	44.5 - 50.8 %
Purity (HPLC)	≥ 92.5 %
¹ H NMR Spectrum	Conforms to Structure

Specification: PRD.0.ZQ5.10000095418

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.