

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
XPhos Pd G4 - 95%

Product Number: **804274**
CAS Number: 1599466-81-5

Formula: C₄₇H₆₄NO₃PPdS
Formula Weight: 860.49 g/mol



TEST

Specification

Appearance (Color)

White to Dark Yellow to Dark Brown and
Grey to Dark Grey

Appearance (Form)

Powder or Crystals

1H NMR Spectrum

Conforms to Structure

31P NMR Spectrum

Conforms to Structure

Carbon Content

62.0 - 69.2 %

Specification: PRD.4.ZQ5.10000077906

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.