

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Acetylcholinesterase from human erythrocytes - buffered aqueous solution, ≥500 units/mg protein (BCA)

Product Number:

C0663

CAS Number:

9000-81-1

Storage Temperature:

2 - 8 °C

TEST

Specification

Appearance (Form)	Conforms
Clear Solution	
Appearance (Colour)	Colorless
Protein Content (Method)	BCA
Enzym. Assay (Spec. Activity)	≥ 500.00 U/mgP
Enzym. Assay (Contam. Name)	Butyrylcholinesterase
Enzym. Assay (Activity U/mgP)	≤ 20.00 U/mgP
Vial Content (Size)	50 - 50
Vial Content (Unit Of Measure)	Units
Vial Content (Actual Content)	50.00 - 75.00
Human Source (HbSAg)	Negative
Human Source (HCV)	Negative
Human Source (HIV)	Negative
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000009958

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.