

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

HAT Media Supplement (50x) Hybri-Max™ -  $\gamma$ -irradiated, lyophilized powder, BioXtra, suitable for hybridoma

Product Number:

H0262

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
10 mL/vial, H <sub>2</sub> O	
Endotoxin Level	≤ 100 EU/ml
when reconstituted at 10 mL/vial	
Sterilized (Gamma Irradiated)	Pass
(Irradiation Dose 25.0 to 35.0 kGy)	
Sterility	Pass
Hybridoma Test	Pass

Specification: PRD.1.ZQ5.10000019859

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.