

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4,4,5,5,6,6,7,7,8,8,9,9,10,10,11,11,11-Heptafluoroundecylamine - ≥97%

Product Number:

4636

CAS Number:

139175-50-1

MDL:


MFCD03788349

Formula:

C₁₁H₈F₁₇N

Formula Weight:

477.16 g/mol

**TEST****Specification**

Appearance (Clarity)

Clear to Turbid

Appearance (Color)

Conforms to Requirements

Colorless/White

Appearance (Form)

Liquid

May Crystallize During Storage

Purity (GC)

≥ 97.0 %

Proton NMR spectrum

Conforms to Structure

Product of Supplier

Confirmed

Fluorous Technologies, Inc. Use of this compound may be protected by U.S. patents 6,156,896; 5,859,247; 5,777,121 and 6,673,539

Specification: PRD.0.ZQ5.10000021791

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.