


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Bromo-5-methylaniline - 97%

Product Number: 686409
CAS Number: 53078-85-6
MDL: MFCD01806278
Formula: C₇H₈BrN
Formula Weight: 186.05 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Yellow to Orange
Appearance (Form)	Liquid or Solid
Infrared spectrum	Conforms to Structure
Purity (Titration by HCL04)	96.5 - 103.5 %
Purity (HPLC)	≥ 96.5 %

Specification: PRD.1.ZQ5.10000007827

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.