

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Platinum on carbon extent of labeling: 5 wt. % loading, matrix activated carbon support

Product Number:

205931**Pt**

MDL:

MFCD00011179

Formula:

Pt

TEST	Specification
Appearance (Color)	Black
Appearance (Form)	Pow der
Loss on Drying Typically 1% to 4%	
Platinum (Pt)	4.0- 5.2 %
ICP Major Analysis Confirms Pt Component	Confirmed

Remarks:

Specification Date : 11/13/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.