


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chlorosulfonylacetyl chloride - 95%

Product Number: 259772
CAS Number: 4025-77-8
MDL: MFCD00007457
Formula: C₂H₂Cl₂O₃S
Formula Weight: 177.01 g/mol



TEST	Specification
Appearance (Color)	Yellow to Brown
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Titration with AgNO ₃	94.5 - 105.5 %

Specification: PRD.0.ZQ5.10000082183

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.