

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Sodium hypochlorite solution - reagent grade, available chlorine 10-15 %

**Product Number:** 425044

NaOCl

MDL: MFCD00011120

Formula: ClNaO

Formula Weight: 74.44 g/mol

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)

Yellow

Appearance (Form)

Liquid

Titration with Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>

10 - 15 %

Available Chlorine by Titration with Sodium Thiosulfate

Recommended Retest Period

1 Year

Specification: PRD.1.ZQ5.10000024089

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.