


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,2-Bis[(2S,5S)-2,5-diethylphospholano]benzene(1,5-cyclooctadiene)rhodium(I) trifluoromethanesulfonate

Product Number: **698407**
CAS Number: 142184-30-3
MDL: MFCD00269861
Formula: C₃₁H₄₈F₃O₃P₂RhS
Formula Weight: 722.62 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Orange to Brown	
Appearance (Form)	Powder
Proton NMR spectrum	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Carbon	49.5 - 53.6 %
Sulfur	4.2 - 4.7 %

Specification: PRD.1.ZQ5.10000009642

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.