


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bradykinin Fragment 1-7 - $\geq 97\%$ (HPLC)

Product Number: B1651
CAS Number: 23815-87-4
MDL: MFCD00076262
Formula: C₃₅H₅₂N₁₀O₉
Formula Weight: 756.85 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 mg/mL, H ₂ O	
Peptide Content	$\geq 70\%$
Purity (HPLC)	$\geq 97\%$

Specification: PRD.0.ZQ5.10000043214

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.