

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl isonipecotate – technical grade

Product Number:

463485

CAS Number:

2971-79-1

MDL:


MFCD00190578

Formula:

C₇H₁₃NO₂

Formula Weight:

143.18 g/mol



TEST

Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Precipitate May Form Upon Standing Which Should Not Affect Application

Proton NMR spectrum

Conforms to Structure

Purity (Titration by HCL04)

≥ 70.0 %

Purity (GC)

≥ 70.0 %

Specification: PRD.1.ZQ5.10000025177

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.