

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Anti-LAMP1–Cy3™ antibody produced in rabbit – 1–2 mg/mL, affinity isolated antibody, buffered aqueous solution
Product Number: L0419

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Purple and Pink
IgG Concentration (OD 280nm)	1.0 - 2.0 mg/ml
F/P Molar Ratio	3.0 - 9.0
ICC/IF Cells	Hela Cells
ICC/IF-Working Concentration	1.0 - 2.0 µg/ml
ICC/IF Cells-2	NRK rat kidney cells
ICC/IF-Working Concentration	1.0 - 2.0 µg/ml
ICC/IF Cells-3	NIH-3T3 Mouse Fibroblasts
ICC/IF-Working Concentration	1.0 - 2.0 µg/ml
ICC/IF Fixation Reagent 1	3%-4% Paraformaldehyde
ICC/IF-Fixation Reagent 2	0.1% Saponin in PBS
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000016928

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.