


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Bromo-1-methylindazole - 97%

Product Number: 717215
CAS Number: 326474-67-3
Formula: C₈H₇BrN₂
Formula Weight: 211.06 g/mol



| TEST | Specification |
|------------------------|--------------------------|
| Appearance (Color) | Conforms to Requirements |
| White to Yellow to Tan | |
| Appearance (Form) | Powder or Solid |
| Infrared spectrum | Conforms to Structure |
| Carbon Content | 43.9 - 47.1 % |
| Nitrogen Content | 12.8 - 13.7 % |
| Purity (GC) | ≥ 96.5 % |

Specification: PRD.1.ZQ5.10000034600

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.