


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Gallamine triethiodide - ≥98% (TLC), powder, muscarinic receptor antagonist

Product Number: G8134
CAS Number: 65-29-2
MDL: MFCD00011832
Formula: C₃₀H₆₀I₃N₃O₃
Formula Weight: 891.53 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Light Brown
Appearance (Form)	Powder
Proton NMR spectrum	Conforms to Structure
Solubility (Color)	Colorless to Beige
Solubility (Turbidity)	Clear to Slightly Hazy
100 mg/mL, H ₂ O	
Carbon	39.4 - 41.2 %
Nitrogen	4.4 - 5.0 %
Purity (TLC)	≥ 98 %

Specification: PRD.1.ZQ5.10000040513

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.