


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,5-Dihydroxynaphthalene - 97%

Product Number: D115606
CAS Number: 83-56-7
MDL: MFCD00003980
Formula: C₁₀H₈O₂
Formula Weight: 160.17 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Purity (HPLC)

Faint Brown to Brown and Grey-Green
Powder
Conforms to Structure
≥ 96.5 %

Specification Date : 04/11/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.