

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Vanillin-(methoxy-d3) - 99 atom % D, 99% (CP)

Product Number:

606170

CAS Number:


74495-74-2

Formula:

C₈D₃H₅O₃

Formula Weight:

155.42 g/mol



TEST

Specification

Appearance (Color)

White to Faint Beige

Appearance (Form)

Powder

Water (by Karl Fischer)

≤ 0.50 %

Purity (GC)

≥ 99.0 %

1H NMR Spectrum

Conforms to Structure

Deuterium NMR Spectrum

Conforms to Structure

D Enrichment

≥ 99.0 %

atom%

Specification: PRD.0.ZQ5.10000102319

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.