

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tellurium dioxide - ≥99%

Product Number: 243450
CAS Number: 7446-07-3
MDL: MFCD00011263
Formula: O₂Te
Formula Weight: 159.6 g/mol



TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

White to Pale Yellow or Green

Appearance (Form)

Powder

X-Ray Diffraction

Conforms to Structure

Gravimetric Analysis

79.1 - 80.8 %

% Te After Reduction With Sulfite

Specification: PRD.1.ZQ5.10000018211

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.