

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Glucuronidase from bovine liver - Type B-3, $\geq 2,000,000$ units/g solid

Product Number:

G0376

CAS Number:

9001-45-0

MDL:

MFCD00130629

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Solubility (Color)

Colorless to Light Yellow

Solubility (Turbidity)

Clear to Slightly Hazy

1 mg/mL in 0.2% NaCl

units/g solid

 ≥ 2000000

One Sigma or modified "Fishman" unit will liberate 1.0 microgram of phenolphthalein from phenolphthalein glucuronide per hour at pH 5.0 at 37 deg C.

Specification: PRD.0.ZQ5.10000036935

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.