


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Di-tert-butylsilyl bis(trifluoromethanesulfonate) - 97%

Product Number: 262021
CAS Number: 85272-31-7
MDL: MFCD00010581
Formula: C₁₀H₁₈F₆O₆S₂Si
Formula Weight: 440.45 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Yellow and Colorless to Light Brown-Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
¹ H NMR Spectrum	Conforms to Structure
Titration with NaOH	96.5 - 103.5 %

Specification: PRD.2.ZQ5.10000031010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.