


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Triphenylphosphine dibromide - 96%

Product Number: **270946**
CAS Number: 1034-39-5
MDL: MFCD00000054
Formula: $C_{18}H_{15}Br_2P$
Formula Weight: 422.09 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Titration by $AgNO_3$
After Oxygen Combustion

White to Light Yellow
Powder or Chunks
Conforms to Structure
95.5 - 104.5 %

Specification: PRD.0.ZQ5.10000051292

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.