

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N,N'-Methylenebisacrylamide – suitable for electrophoresis (after filtration or allowing insolubles to settle)

Product Number:

M7256

CAS Number:

110-26-9

MDL:


MFCD00008625

Formula:

C₇H₁₀N₂O₂

Formula Weight:

154.17 g/mol



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Solubility (Color)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear to Turbid

80 mg/ml, MeOH

Infrared spectrum

Conforms to Structure

Carbon

53.4 - 55.6 %

Nitrogen

17.8 - 18.5 %

Suitability

Suitable

Suitable for use in electrophoresis after filtration or allowing insolubles to settle

Specification: PRD.0.ZQ5.10000036288

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.