

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(S)[(Sp)-2-(Diphenylphosphino)ferrocenyl]-4-isopropylloxazoline - 97%

Product Number: **717398**
CAS Number: 163169-29-7

Formula: C₂₈H₂₈FeNOP
Formula Weight: 481.35 g/mol



TEST	Specification
Appearance (Color)	Yellow to Orange to Brown
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
¹ H NMR Spectrum	Conforms to Structure
Carbon Content	67.40 - 72.30 %
Optical Rotation	15 - 48 deg
in chloroform w/amylene	
Product of Supplier	Confirmed
NAMENA CORP	

Specification: PRD.0.ZQ5.10000101954

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.