


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(R)-DM-BINAP

Product Number: **692379**
CAS Number: 137219-86-4
MDL: MFCD01630821
Formula: C₅₂H₄₈P₂
Formula Weight: 734.89 g/mol



TEST

Specification

Appearance (Color)	White to Beige
Appearance (Form)	Solid or Powder or Crystals or Chunks
Proton NMR spectrum	Conforms to Structure
Proton NMR	≥ 94.0 %
Optical Rotation	164.0 - 174.0 deg
c = 1%; Chloroform	
Product of Supplier	Confirmed
Takasago International Corporation	

Specification: PRD.1.ZQ5.10000037400

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.