


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium iodate - $\geq 99\%$

Product Number: **S4007**
CAS Number: 7681-55-2
MDL: MFCD00003533
Formula: INaO_3
Formula Weight: 197.89 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Turbidity) c = 50mg/ml; Water	Clear to Very Slightly Hazy
Solubility (Color)	Colorless
Titration by $\text{Na}_2\text{S}_2\text{O}_3$	$\geq 99.0\%$

Specification: PRD.0.ZQ5.10000030496

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.