


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Isethionic acid sodium salt - 98%

Product Number: 220078
CAS Number: 1562-00-1
MDL: MFCD00007534
Formula: C₂H₅NaO₄S
Formula Weight: 148.11 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Conforms to Requirements
Powder, Crystalline Powder or Granular Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	97.5 - 102.5 %
After Ion Exchange	
ICP Major Analysis	Confirmed
Confirms Sodium Component	
Solubility (Turbidity)	Clear
C = 10% in H ₂ O	
Solubility (Color)	Colorless

Specification: PRD.1.ZQ5.10000008922

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.