


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4,6-Dichlororesorcinol - 97%

Product Number: 130982
CAS Number: 137-19-9
MDL: MFCD00002274
Formula: C₆H₄Cl₂O₂
Formula Weight: 179.00 g/mol



TEST

Specification

Appearance (Color)

White to Grey and Faint Beige to Beige

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000042337

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.