


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Fast Red Violet LB base

Product Number: **274054**
CAS Number: 121-22-2
MDL: MFCD00007780
Formula: C₁₄H₁₃ClN₂O
Formula Weight: 260.72 g/mol



TEST

Specification

Appearance (Color)	White to Light Yellow to Light Brown
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	62.0 - 67.0 %
Nitrogen	10.0 - 12.0 %
Solubility (Turbidity)	Clear
2.5% in Acetone	
Solubility (Color)	Colorless to Yellow to Brown
2.5% in Acetone	

Specification: PRD.0.ZQ5.10000054362

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.