

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Cycloheximide solution – Ready-Made Solution, microbial, 100 mg/mL in DMSO, 0.2 µm filtered

Product Number:

C4859

Formula:


C<sub>15</sub>H<sub>23</sub>NO<sub>4</sub>

Formula Weight:

281.35 g/mol

Storage Temperature:

2 - 8 °C



### TEST

### Specification

Appearance (Form)

Clear Solution

Conforms

Appearance (Colour)

Concentration

Suitability (Result)

Suitability (Test)

Recommended Retest Period

2 Years

Light Yellow to Yellow

100.000 - 105.000 mg/ml

Pass

CHO C. Growth Inhib.

-----

Specification: PRE.0.ZQ5.10000009960

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.