

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Copper(I) iodide - 98%

Product Number: 205540
CAS Number: 7681-65-4

Formula: CuI
Formula Weight: 190.45 g/mol

CuI

TEST	Specification
------	---------------

Appearance (Colour)	White to Grey and Faint Beige to Beige and Faint Brown to Brown and Grey-Brown
Appearance (Form)	Powder
Redox Titration	97.5 - 102.5 %
ICP Major Analysis	Corresponds to Requirements
Copper Component: Corresponds	

Specification: PRE.0.ZQ5.10000018163

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.