

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tributyl phosphate - 97%

Product Number: 158615
CAS Number: 126-73-8
MDL: MFCD00009436
Formula: C₁₂H₂₇O₄P
Formula Weight: 266.31 g/mol



TEST

Specification

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 20 ° C	1.422 - 1.426
Purity (GC)	≥ 96.5 %
Recommended Retest Period	-----
2.5 Years	

Specification Date : 04/20/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.