

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
R-Alpine-Hydride® - 0.5 M in THF

Product Number: 229024

MDL: MFCD00082426

Formula: C₁₈H₃₂BLi

TEST	Specification
------	---------------

Appearance (Color)	Colorless to Very Faint Yellow
Appearance (Form)	Liquid
Optical Rotation (Degrees) NEAT	-10.0 - -7.0
Hydrogen Evolution	0.40 - 0.60 M

Specification: PRD.0.ZQ5.10000021908

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.