

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Flumequine - ≥97.0% (HPLC)

Product Number: **F7016**
CAS Number: 42835-25-6

Formula: C₁₄H₁₂FNO₃
Formula Weight: 261.3 AMU



TEST	Specification
Appearance (Form) Powder	Conforms
Appearance (Color) White to Off-White	Conforms
Purity (HPLC)	≥ 97.0 %
Identification (¹ H-NMR)	Conforms
Heavy Metals < = 10ppm	Conforms
Loss on Drying	≤ 1.0 %
Sulphated Ash	≤ 0.10 %
Recommended Retest Period 2Years	-----

Specification: PRE.1.ZQ5.10000014690

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.