

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium percarbonate - avail. H₂O₂ 20-30 %**Product Number:****371432**Na₂CO₃ · 1.5 H₂O₂

CAS Number:

15630-89-4

MDL:

MFCD00043204

Formula:

CNa₂O₃ · 1.5H₂O₂

Formula Weight:

157.01 g/mol

TEST**Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Conforms

Powder, Beads or Granular Powder

ICP Major Analysis:

Confirmed

Confirms Sodium Component

Titration by HCL (% CO₃)Typically Ca. 35% - 38% CO₃Titration by Na₂S₂O₃Approximately 20% - 30% as H₂O₂

pH

10.0 - 11.0

c = 3%; Water

Specification: PRD.0.ZQ5.10000043723

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.