

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Titanium - foil, thickness 0.127 mm, ≥99.99% trace metals basis

Product Number: 460397
CAS Number: 7440-32-6
MDL: MFCD00011264
Formula: Ti
Formula Weight: 47.87 g/mol

Ti

TEST	Specification
Appearance (Color)	Silver
Appearance (Form)	Foil
X-Ray Diffraction	Conforms to Structure
Thickness	Conforms
0.127 mm	
Trace Metal Analysis	≤ 100.0 ppm
Purity	Meets Requirements
> = 99.99% Based On Trace Metals Analysis	

Specification: PRD.0.ZQ5.10000038532

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.