


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Cholecalciferol - meets USP testing specifications

**Product Number:** C1357  
**CAS Number:** 67-97-0  
**MDL:** MFCD00078131  
**Formula:** C<sub>27</sub>H<sub>44</sub>O  
**Formula Weight:** 384.64 g/mol  
**Storage Temperature:** 2 - 8 °C



TEST	Specification
Identity	Pass
Specific Rotation	105 - 112 °
Assay	97.0 - 103.0 %
Residual Solvents	Meets Requirements

Specification: PRD.2.ZQ5.10000055777

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.