

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anti-Protein A antibody produced in rabbit - fractionated antiserum, lyophilized powder

Product Number: P3775

TEST	Specification
Appearance (Form)	Powder
Appearance (Color)	White to Off White
Solubility (Solvent)	Water
Solubility (Volume)	2 - 2 ml
Solubility (Color)	Colorless to Slightly Yellow
Solubility (Turbidity)	Clear
Protein by Biuret	4 - 100 mg/ml
ELISA - Titer	1:100,000-1:200,000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000019691

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.