


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Product Number: **795127**Formula: C<sub>10</sub>H<sub>13</sub>BF<sub>3</sub>KN<sub>2</sub>O<sub>2</sub>

Formula Weight: 300.13 g/mol

**TEST****Specification**

Appearance (Color)	White to Faint Yellow and Faint Beige and Faint Brown
Appearance (Form)	Powder or Crystals
<sup>1</sup> H NMR Spectrum	Conforms to Structure
<sup>19</sup> F NMR Spectrum	Conforms to Structure
<sup>11</sup> B NMR Spectrum	Conforms to Structure
Carbon	
Specification: Report Result	
Nitrogen	
Specification: Report Result	
Purity (HPLC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000061568

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.